

AeroNOx 2.0™

Nitric Oxide Delivery

Nitric Oxide Delivery System



NICU



Air



Ground

Meets transport standards
Lightweight and rugged
NO, NO₂ and O₂ monitoring

AeroNOx 2.0™

Nitric Oxide Delivery

Standard Features



LCD Color Screen

Easy to read and interpret color display.



Easy Maintenance

Easy access to replace and service the battery and sensors.



Smart Alarms

Visual and audible alarms updated to the latest alarm standards.

Transport Mount



Order Information

Part #	Component	Part #	Component
731-0317	AeroNOx 2.0 Nitric Oxide Delivery System	731-9147	INostat Kit
738-1861	AeroNOx 2.0 Delivery Hose, 3ft, Stainless Steel	731-0274	AeroNOx 2.0 Calibration Kit
888-0115	AeroNOx 2.0 Rechargeable Battery	731-0276	AeroNOx 2.0 Sample/Delivery Kit, Conventional (5/pack)
731-0272	AeroNOx 2.0 Transport Mount	731-0277	AeroNOx 2.0 Sample/Delivery Kit, HFV (5/pack)
731-0350	AeroNOx 2.0 Pole Clamp Assembly	731-0278	AeroNOx 2.0 Bagger (5/pack)
731-9919	INostat Baggers (5/pack)	738-1915	4' Power Cord (AeroNOx 2.0 to Incubator)
317-0735	Clamp and Cradle Kit (2 pair/kit)	738-1914	1.5' Power Cord (AeroNOx 2.0 to Incubator)

Europe

International Biomedical AG
Glütschbachstrasse 2, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical
8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.